

MHC class 1 Conjugated Antibody

Catalog No: #C56034



Package Size: #C56034-AF350 100ul #C56034-AF405 100ul #C56034-AF488 100ul #C56034-AF555 100ul #C56034-AF594 100ul #C56034-AF647 100ul #C56034-AF680 100ul #C56034-AF750 100ul #C56034-Biotin 100ul #C56034-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	MHC class 1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	MHC class 1 Antibody detects endogenous levels of total MHC class 1
Immunogen Description	A synthesized peptide derived from human MHC class 1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Aw-68;HLA class I histocompatibility antigen, A-28 alpha chain;MHC class I antigen A*68;HLA-A; MHC class I antigen HLA A heavy chain;
Accession No.	Uniprot:P01891
Calculated MW	41kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

MHC class 1 is involved in the presentation of foreign antigens to the immune system.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.