

CD47 Conjugated Antibody

Catalog No: #C56041



Package Size: #C56041-AF350 100ul #C56041-AF405 100ul #C56041-AF488 100ul #C56041-AF555 100ul #C56041-AF594 100ul #C56041-AF647 100ul #C56041-AF680 100ul #C56041-AF750 100ul #C56041-Biotin 100ul #C56041-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CD47 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	CD47 Antibody detects endogenous levels of total CD47
Immunogen Description	A synthesized peptide derived from human CD47
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CD47; IAP; Integrin Associated Protein; MER6 ; OA3; Protein MER6 ; Antigen identified by monoclonal antibody 1D8
Accession No.	Uniprot:Q08722
Calculated MW	55kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

FC: 1:50-1:200

Product Description

CD47 is an integral membrane protein that plays a role in the regulation of cation fluxes across cell membranes. Specifically, CD47 is involved in the increase in intracellular calcium concentration that occurs upon cell adhesion to the extracellular matrix. It is also a receptor for the C-terminal cell binding domain of thrombospondin (SIRP). CD47 is absent from Rh-null erythrocytes, but does play a role in cell adhesion in non-erythroid cells and may prevent premature elimination of erythrocytes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.