

PELP1 Conjugated Antibody

Catalog No: #C56045

Package Size: #C56045-AF350 100ul #C56045-AF405 100ul #C56045-AF488 100ul #C56045-AF555 100ul #C56045-AF594 100ul #C56045-AF647 100ul #C56045-AF680 100ul #C56045-AF750 100ul #C56045-Biotin 100ul #C56045-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PELP1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	PELP1 Antibody detects endogenous levels of total PELP1
Immunogen Description	A synthesized peptide derived from human PELP1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	HMX3; MNAR; P160; PELP1; PELP1 proline glutamic acid leucine rich protein 1; PELP1 proline- glutamic;
Accession No.	Uniprot:Q8IZL8
Calculated MW	160kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

This gene encodes a transcription factor which coactivates transcription of estrogen receptor responsive genes and corepresses genes activated by other hormone receptors or sequence-specific transcription factors. Expression of this gene is regulated by both members of the estrogen receptor family. This gene may be involved in the progression of several types of cancer. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.