

alpha + beta Synuclein Conjugated Antibody

Catalog No: #C56061

Package Size: #C56061-AF350 100ul #C56061-AF405 100ul #C56061-AF488 100ul #C56061-AF555 100ul #C56061-AF594 100ul #C56061-AF647 100ul #C56061-AF680 100ul #C56061-AF750 100ul #C56061-Biotin 100ul #C56061-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | alpha + beta Synuclein Conjugated Antibody |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB, IF, FC |
| Species Reactivity | Human Mouse Rat |
| Specificity | alpha + beta Synuclein Antibody detects endogenous levels of total alpha + beta Synuclein |
| Immunogen Description | A synthesized peptide derived from human alpha + beta Synuclein |
| Conjugates | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750 |
| Other Names | alpha synuclein; beta synuclein ; NACP; PD1 ; SNCA; SNCB; |
| Accession No. | Uniprot:P37840/Q16143 |
| Calculated MW | 18kDa |
| Storage | Store at 4°C in dark for 6 months |

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Alpha-synuclein is a member of the synuclein family, which also includes beta- and gamma-synuclein. Synucleins are abundantly expressed in the brain and alpha- and beta-synuclein inhibit phospholipase D2 selectively.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.