

# beta II Tubulin Conjugated Antibody

Catalog No: #C56067

Package Size: #C56067-AF350 100ul #C56067-AF405 100ul #C56067-AF488 100ul #C56067-AF555 100ul #C56067-AF594 100ul #C56067-AF647 100ul #C56067-AF680 100ul #C56067-AF750 100ul #C56067-Biotin 100ul #C56067-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	beta II Tubulin Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	beta II Tubulin Antibody detects endogenous levels of total beta II Tubulin
Immunogen Description	A synthesized peptide derived from human beta II Tubulin
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	class II beta tubulin isotype; Beta2; Tubulin beta 2 chain; TUBB; TUBB2;
Accession No.	Uniprot:P68371/Q13885/Q9BVA1
Calculated MW	50kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.