

# Transferrin Conjugated Antibody

Catalog No: #C56069

Package Size: #C56069-AF350 100ul #C56069-AF405 100ul #C56069-AF488 100ul #C56069-AF555 100ul #C56069-AF594 100ul #C56069-AF647 100ul #C56069-AF680 100ul #C56069-AF750 100ul #C56069-Biotin 100ul #C56069-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Transferrin Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	Transferrin Antibody detects endogenous levels of total Transferrin
Immunogen Description	A synthesized peptide derived from human Transferrin
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	TFQTL1; TF; TRFE; Transferrin; Serotransferrin; Serotransferrin;
Accession No.	Uniprot:P02787
Calculated MW	77kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Transferrins are iron binding transport proteins which can bind two Fe<sup>3+</sup> ions in association with the binding of an anion, usually bicarbonate. It is responsible for the transport of iron from sites of absorption and heme degradation to those of storage and utilization. Serum transferrin may also have a further role in stimulating cell proliferation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.