

Brd4 Conjugated Antibody

Catalog No: #C56072

Package Size: #C56072-AF350 100ul #C56072-AF405 100ul #C56072-AF488 100ul #C56072-AF555 100ul #C56072-AF594 100ul #C56072-AF647 100ul #C56072-AF680 100ul #C56072-AF750 100ul #C56072-Biotin 100ul #C56072-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Brd4 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Brd4 Antibody detects endogenous levels of total Brd4
Immunogen Description	A synthesized peptide derived from human Brd4
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Brd4; CAP; HUNK1; MCAP; Bromodomain containing 4; chromosome associated protein;
Accession No.	Uniprot:O60885
Calculated MW	152kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Bromodomain-containing protein 4 (BRD4) is a member of the bromodomains and extra terminal (BET) family of proteins, which also includes BRD2, BRD3, and BRDT. BET family proteins contain two tandem bromodomains and an extra terminal (ET) domain, and bind acetyl lysine residues. BRD4 is a chromatin-binding protein with a preference for Lys14 on histone H3 as well as Lys5 and Lys12 on histone H4. BRD4 chromatin binding occurs throughout the cell cycle, including condensed mitotic chromosomes, when the majority of genes are silenced.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.