

JNK1/2/3 Conjugated Antibody

Catalog No: #C56078

Package Size: #C56078-AF350 100ul #C56078-AF405 100ul #C56078-AF488 100ul #C56078-AF555 100ul #C56078-AF594 100ul #C56078-AF647 100ul #C56078-AF680 100ul #C56078-AF750 100ul #C56078-Biotin 100ul #C56078-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	JNK1/2/3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat Cow Dog Monkey Zebrafish
Specificity	JNK1+JNK2+JNK3 Antibody detects endogenous levels of total JNK1+JNK2+JNK3
Immunogen Description	A synthesized peptide derived from human JNK1+JNK2+JNK3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	c-Jun N-terminal kinase 3; JNK1; JNK2;JNK3; JNK3 alpha protein kinase; JNK3A; MAP kinase 10; MAP kinase p49 3F12; MAPK 10; MAPK10;
Accession No.	Uniprot:P45983/P45984/P53779
Calculated MW	46,54kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Serine/threonine-protein kinase involved in various processes such as cell proliferation, differentiation, migration, transformation and programmed cell death. Extracellular stimuli such as proinflammatory cytokines or physical stress stimulate the stress-activated protein kinase/c-Jun N-terminal kinase (SAP/JNK) signaling pathway.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.