

GATA2/3 Conjugated Antibody

Catalog No: #C56082



Package Size: #C56082-AF350 100ul #C56082-AF405 100ul #C56082-AF488 100ul #C56082-AF555 100ul #C56082-AF594 100ul
 #C56082-AF647 100ul #C56082-AF680 100ul #C56082-AF750 100ul #C56082-Biotin 100ul #C56082-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GATA2/3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	GATA2/3 Antibody detects endogenous levels of total GATA2/3
Immunogen Description	A synthesized peptide derived from human GATA2/3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	DCML;GATA binding protein 2;GATA binding protein 3;HDR;HDRS;IMD21;MONOMAC;NFE1B
Accession No.	Uniprot:P23769/P23771
Calculated MW	51kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Transcriptional activator which regulates endothelin-1 gene expression in endothelial cells. Binds to the consensus sequence 5'-AGATAG-3'.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.