

# KDM1/LSD1 Conjugated Antibody

Catalog No: #C56086

Package Size: #C56086-AF350 100ul #C56086-AF405 100ul #C56086-AF488 100ul #C56086-AF555 100ul #C56086-AF594 100ul #C56086-AF647 100ul #C56086-AF680 100ul #C56086-AF750 100ul #C56086-Biotin 100ul #C56086-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	KDM1/LSD1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	KDM1 / LSD1 Antibody detects endogenous levels of total KDM1 / LSD1
Immunogen Description	A synthesized peptide derived from human KDM1 / LSD1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	KDM1;AOF2;BHC110;KIAA0601;LSD1; KDM1A;
Accession No.	Uniprot:O60341
Calculated MW	110kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Histone demethylase that demethylates both 'Lys-4' (H3K4me) and 'Lys-9' (H3K9me) of histone H3, thereby acting as a coactivator or a corepressor, depending on the context. Acts by oxidizing the substrate by FAD to generate the corresponding imine that is subsequently hydrolyzed. Acts as a corepressor by mediating demethylation of H3K4me, a specific tag for epigenetic transcriptional activation. Demethylates both mono- (H3K4me1) and di-methylated (H3K4me2) H3K4me. May play a role in the repression of neuronal genes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.