

# Lin28B Conjugated Antibody

Catalog No: #C56087

Package Size: #C56087-AF350 100ul #C56087-AF405 100ul #C56087-AF488 100ul #C56087-AF555 100ul #C56087-AF594 100ul #C56087-AF647 100ul #C56087-AF680 100ul #C56087-AF750 100ul #C56087-Biotin 100ul #C56087-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Lin28B Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	Lin28B Antibody detects endogenous levels of total Lin28B
Immunogen Description	A synthesized peptide derived from human Lin28B
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CSDD2; FLJ16517; LIN28B; Protein lin-28 homolog B;
Accession No.	Uniprot:Q6ZN17
Calculated MW	30kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

It has recently been shown that upregulation of LIN28A and LIN28B in primary human tumors and human cancer cell lines is correlated with downregulation of let-7 miRNAs. LIN28 genes are reported to be involved in primordial germ cell development and germ cell malignancy. In addition, allelic variation in LIN28B is associated with regulating the timing of puberty in humans.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.