

# Hsp90 Conjugated Antibody

Catalog No: #C56090

Package Size: #C56090-AF350 100ul #C56090-AF405 100ul #C56090-AF488 100ul #C56090-AF555 100ul #C56090-AF594 100ul #C56090-AF647 100ul #C56090-AF680 100ul #C56090-AF750 100ul #C56090-Biotin 100ul #C56090-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Hsp90 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	Hsp90 Antibody detects endogenous levels of total Hsp90
Immunogen Description	A synthesized peptide derived from human Hsp90
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	HSP 90; HSP90A; HSP90B; LAP-2; HSPC1; HSPC2;
Accession No.	Uniprot:P07900/P08238
Calculated MW	90kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Molecular chaperone that promotes the maturation, structural maintenance and proper regulation of specific target proteins involved for instance in cell cycle control and signal transduction. Undergoes a functional cycle that is linked to its ATPase activity. This cycle probably induces conformational changes in the client proteins, thereby causing their activation.

**Note:** This product is for in vitro research use only and is not intended for use in humans or animals.