

CCR3 Conjugated Antibody

Catalog No: #C56092

Package Size: #C56092-AF350 100ul #C56092-AF405 100ul #C56092-AF488 100ul #C56092-AF555 100ul #C56092-AF594 100ul #C56092-AF647 100ul #C56092-AF680 100ul #C56092-AF750 100ul #C56092-Biotin 100ul #C56092-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CCR3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	CCR3 Antibody detects endogenous levels of total CCR3
Immunogen Description	A synthesized peptide derived from human CCR3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CCR3; CD193; CCCKR3; CMKBR3; CKR3; CMKBR1L2; C-C chemokine receptor type 3;
Accession No.	Uniprot:P51677
Calculated MW	60kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Receptor for a C-C type chemokine. Binds to eotaxin, eotaxin-3, MCP-3, MCP-4, RANTES and MIP-1 delta. Subsequently transduces a signal by increasing the intracellular calcium ions level. Alternative coreceptor with CD4 for HIV-1 infection.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.