

Cyclin A1/A2 Conjugated Antibody

Catalog No: #C56107

Package Size: #C56107-AF350 100ul #C56107-AF405 100ul #C56107-AF488 100ul #C56107-AF555 100ul #C56107-AF594 100ul #C56107-AF647 100ul #C56107-AF680 100ul #C56107-AF750 100ul #C56107-Biotin 100ul #C56107-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Cyclin A1/A2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Cyclin A1/A2 Antibody detects endogenous levels of Cyclin A1/A2
Immunogen Description	A synthesized peptide derived from human Cyclin A1/A2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CCN1; CCNA; CCNA1; CCNA2; CT146; Cyclin-A; Cyclin-A1; Cyclin-A2;
Accession No.	Uniprot:P20248/P78396
Calculated MW	52,49kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

May primarily function in the control of the germline meiotic cell cycle and additionally in the control of mitotic cell cycle in some somatic cells. Essential for the control of the cell cycle at the G1/S (start) and the G2/M (mitosis) transitions.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.