

STAT5A/B Conjugated Antibody

Catalog No: #C56109



Package Size: #C56109-AF350 100ul #C56109-AF405 100ul #C56109-AF488 100ul #C56109-AF555 100ul #C56109-AF594 100ul #C56109-AF647 100ul #C56109-AF680 100ul #C56109-AF750 100ul #C56109-Biotin 100ul #C56109-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	STAT5A/B Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	STAT5A/B Antibody detects endogenous levels of STAT5A/B
Immunogen Description	A synthesized peptide derived from human STAT5A/B
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Mammary gland factor (STAT5A); Signal transducer and activator of transcription 5A; Stat 5a; STAT 5b; STAT5A; Transcription factor STAT5B;
Accession No.	Uniprot:P42229/P51692
Calculated MW	90kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

FC: 1:50-1:200

Product Description

Carries out a dual function: signal transduction and activation of transcription. Binds to the GAS element and activates PRL-induced transcription.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.