

MEK3/MEK6 Conjugated Antibody

Catalog No: #C56110

Package Size: #C56110-AF350 100ul #C56110-AF405 100ul #C56110-AF488 100ul #C56110-AF555 100ul #C56110-AF594 100ul
#C56110-AF647 100ul #C56110-AF680 100ul #C56110-AF750 100ul #C56110-Biotin 100ul #C56110-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MEK3/MEK6 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	MEK3/MEK6 Antibody detects endogenous levels of MEK3/MEK6
Immunogen Description	A synthesized peptide derived from human MEK3/MEK6
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	MEK6; MKK6; MAPKK6; PRKMK6; SAPKK3; MAP2K6; MEK3; MAP kinase kinase 3; MAPKK3; MAPK/ERK kinase 3;
Accession No.	Uniprot:P46734/P52564
Calculated MW	39,37kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

FC: 1:50-1:200

Product Description

Dual specificity kinase. Is activated by cytokines and environmental stress in vivo. Catalyzes the concomitant phosphorylation of a threonine and a tyrosine residue in the MAP kinase p38.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.