

CDCA3 Conjugated Antibody

Catalog No: #C56115

Package Size: #C56115-AF350 100ul #C56115-AF405 100ul #C56115-AF488 100ul #C56115-AF555 100ul #C56115-AF594 100ul #C56115-AF647 100ul #C56115-AF680 100ul #C56115-AF750 100ul #C56115-Biotin 100ul #C56115-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CDCA3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	CDCA3 Antibody detects endogenous levels of total CDCA3
Immunogen Description	A synthesized peptide derived from human CDCA3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Cell division cycle-associated protein 3; Trigger of mitotic entry protein 1; TOME-1; Gene-rich cluster protein C8; CDCA3; C8; GRCC8; TOME1;
Accession No.	Uniprot:Q99618
Calculated MW	26kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

F-box-like protein which is required for entry into mitosis. Acts by participating in E3 ligase complexes that mediate the ubiquitination and degradation of WEE1 kinase at G2/M phase (By similarity).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.