

AKT1/2 Conjugated Antibody

Catalog No: #C56120

Package Size: #C56120-AF350 100ul #C56120-AF405 100ul #C56120-AF488 100ul #C56120-AF555 100ul #C56120-AF594 100ul #C56120-AF647 100ul #C56120-AF680 100ul #C56120-AF750 100ul #C56120-Biotin 100ul #C56120-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	AKT1/2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	AKT1/2 Antibody detects endogenous levels of total AKT1/2
Immunogen Description	A synthesized peptide derived from human AKT1/2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AKT1; AKT1 kinase; AKT2; AKT2 kinase; C AKT; CWS6; Oncogene AKT1; PKB alpha; PKB beta; PRKBA; PRKBB;
Accession No.	Uniprot:P31749/P31751
Calculated MW	56kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

FC: 1:50-1:200

Product Description

The serine/threonine kinase AKT (protein kinase B or PKB) has a central role in the regulation of several signaling pathways controlling cell proliferation, apoptosis, angiogenesis, and diabetes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.