

iNOS Conjugated Antibody

Catalog No: #C56122

Package Size: #C56122-AF350 100ul #C56122-AF405 100ul #C56122-AF488 100ul #C56122-AF555 100ul #C56122-AF594 100ul #C56122-AF647 100ul #C56122-AF680 100ul #C56122-AF750 100ul #C56122-Biotin 100ul #C56122-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	iNOS Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse
Specificity	iNOS Antibody detects endogenous levels of total iNOS
Immunogen Description	A synthesized peptide derived from human iNOS
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	HEP-NOS; Hepatocyte NOS; Inducible NO synthase; Inducible NOS; iNOS; NOS2;
Accession No.	Uniprot:P35228
Calculated MW	131kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

iNOS Produces nitric oxide (NO) which is a messenger molecule with diverse functions throughout the body. In macrophages, NO mediates tumoricidal and bactericidal actions.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.