

CD3 zeta Conjugated Antibody

Catalog No: #C56124

Package Size: #C56124-AF350 100ul #C56124-AF405 100ul #C56124-AF488 100ul #C56124-AF555 100ul #C56124-AF594 100ul #C56124-AF647 100ul #C56124-AF680 100ul #C56124-AF750 100ul #C56124-Biotin 100ul #C56124-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CD3 zeta Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	IF, FC
Species Reactivity	Human
Specificity	CD3 zeta Antibody detects endogenous levels of total CD3 zeta
Immunogen Description	A synthesized peptide derived from human CD3 zeta
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CD3-zeta; T-cell receptor T3 zeta chain; T3Z; TCRZ; CD247;
Accession No.	Uniprot:P20963
Calculated MW	kDa
Storage	Store at 4°C in dark for 6 months

Application Details

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Defects in CD3D are a cause of severe combined immunodeficiency autosomal recessive T-cell-negative/B-cell-positive/NK-cell-positive (T(-)/B(+)/NK(+)) SCID [MIM:608971]. A form of severe combined immunodeficiency (SCID), a genetically and clinically heterogeneous group of rare congenital disorders characterized by impairment of both humoral and cell-mediated immunity, leukopenia, and low or absent antibody levels.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.