

BMP11 Conjugated Antibody

Catalog No: #C56128

Package Size: #C56128-AF350 100ul #C56128-AF405 100ul #C56128-AF488 100ul #C56128-AF555 100ul #C56128-AF594 100ul #C56128-AF647 100ul #C56128-AF680 100ul #C56128-AF750 100ul #C56128-Biotin 100ul #C56128-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	BMP11 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Rat
Specificity	BMP11 Antibody detects endogenous levels of total BMP11
Immunogen Description	A synthesized peptide derived from human BMP11
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	GDF11; BMP11; Bone morphogenetic protein 11; GDF-11; BMP-11;
Accession No.	Uniprot:O95390
Calculated MW	45kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Secreted signal that acts globally to specify positional identity along the anterior/posterior axis during development. Play critical roles in patterning both mesodermal and neural tissues and in establishing the skeletal pattern.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.