

Myosin heavy chain Conjugated Antibody

Catalog No: #C56130

Package Size: #C56130-AF350 100ul #C56130-AF405 100ul #C56130-AF488 100ul #C56130-AF555 100ul #C56130-AF594 100ul
#C56130-AF647 100ul #C56130-AF680 100ul #C56130-AF750 100ul #C56130-Biotin 100ul #C56130-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myosin heavy chain Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Myosin heavy chain Antibody detects endogenous levels of total Myosin heavy chain
Immunogen Description	A synthesized peptide derived from human Myosin heavy chain
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	MYH6; MYH7; MYHC; MYHCA; MYHCB; Myosin heavy chain 6; Myosin heavy chain; Myosin-6;
Accession No.	Uniprot:P12883/P13533
Calculated MW	224kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Muscle contraction.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.