

## PDI Conjugated Antibody

Catalog No: #C56131



Package Size: #C56131-AF350 100ul #C56131-AF405 100ul #C56131-AF488 100ul #C56131-AF555 100ul #C56131-AF594 100ul #C56131-AF647 100ul #C56131-AF680 100ul #C56131-AF750 100ul #C56131-Biotin 100ul #C56131-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PDI Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	PDI Antibody detects endogenous levels of total PDI
Immunogen Description	A synthesized peptide derived from human PDI
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	P4HB; PDIA2; PDIR;
Accession No.	Uniprot:Q13087
Calculated MW	58kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Acts as an intracellular estrogen-binding protein. May be involved in modulating cellular levels and biological functions of estrogens in the pancreas. May act as a chaperone that inhibits aggregation of misfolded proteins.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.