

Ras Conjugated Antibody

Catalog No: #C56133

Package Size: #C56133-AF350 100ul #C56133-AF405 100ul #C56133-AF488 100ul #C56133-AF555 100ul #C56133-AF594 100ul #C56133-AF647 100ul #C56133-AF680 100ul #C56133-AF750 100ul #C56133-Biotin 100ul #C56133-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Ras Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Ras Antibody detects endogenous levels of total Ras
Immunogen Description	A synthesized peptide derived from human Ras
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	GTPase HRas; GTPase KRas; GTPase NRas; H-Ras-1; H-RASIDX; Ha-Ras; HAMSIV; HRAS; K RAS2A; K RAS2B; K RAS4A; K RAS4B; K-RAS; KRAS; N-RAS; N-terminally processed; NRAS; p21ras; RASH1; RASK2;
Accession No.	Uniprot:P01111/P01112/P01116
Calculated MW	22kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

Product Description

The 21 kDa guanine-nucleotide binding proteins (K-Ras, H-Ras, and N-Ras) cycle between active (GTP-bound) and inactive (GDP-bound) forms. Ras proteins bind GDP/GTP and possess intrinsic GTPase activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.