

AQP4 Conjugated Antibody

Catalog No: #C56144



Package Size: #C56144-AF350 100ul #C56144-AF405 100ul #C56144-AF488 100ul #C56144-AF555 100ul #C56144-AF594 100ul

#C56144-AF647 100ul #C56144-AF680 100ul #C56144-AF750 100ul #C56144-Biotin 100ul #C56144-Conjugated 50ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	AQP4 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	IF
Species Reactivity	Human Rat
Specificity	AQP4 Antibody detects endogenous levels of total AQP4
Immunogen Description	A synthesized peptide derived from human AQP4
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AQP4; Aquaporin type 4; HMIWC2; Mercurial insensitive water channel; MIWC; WCH4;
Accession No.	Uniprot:P55087
Calculated MW	kDa
Storage	Store at 4°C in dark for 6 months

Application Details

IF:1:50-1:200

Product Description

Forms a water-specific channel. Osmoreceptor which regulates body water balance and mediates water flow within the central nervous system.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.