

BMP7 Conjugated Antibody

Catalog No: #C56145



Package Size: #C56145-AF350 100ul #C56145-AF405 100ul #C56145-AF488 100ul #C56145-AF555 100ul #C56145-AF594 100ul #C56145-AF647 100ul #C56145-AF680 100ul #C56145-AF750 100ul #C56145-Biotin 100ul #C56145-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	BMP7 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	BMP7 Antibody detects endogenous levels of total BMP7
Immunogen Description	A synthesized peptide derived from human BMP7
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Bmp7; Bone morphogenetic protein 7; Eptotermin alfa; OP 1; Osteogenic protein 1;
Accession No.	Uniprot:P18075
Calculated MW	49kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Induces cartilage and bone formation. May be the osteoinductive factor responsible for the phenomenon of epithelial osteogenesis. Plays a role in calcium regulation and bone homeostasis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.