

FABP1 Conjugated Antibody

Catalog No: #C56146

Package Size: #C56146-AF350 100ul #C56146-AF405 100ul #C56146-AF488 100ul #C56146-AF555 100ul #C56146-AF594 100ul #C56146-AF647 100ul #C56146-AF680 100ul #C56146-AF750 100ul #C56146-Biotin 100ul #C56146-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	FABP1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse
Specificity	FABP1 Antibody detects endogenous levels of total FABP1
Immunogen Description	A synthesized peptide derived from human FABP1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	FABP1; FABPL; Fatty acid binding protein 1 liver; L FABP;
Accession No.	Uniprot:P07148
Calculated MW	14kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Binds free fatty acids and their coenzyme A derivatives, bilirubin, and some other small molecules in the cytoplasm. May be involved in intracellular lipid transport.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.