

Cathepsin L/V/K/H Conjugated Antibody

Catalog No: #C56147

Package Size: #C56147-AF350 100ul #C56147-AF405 100ul #C56147-AF488 100ul #C56147-AF555 100ul #C56147-AF594 100ul #C56147-AF647 100ul #C56147-AF680 100ul #C56147-AF750 100ul #C56147-Biotin 100ul #C56147-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Cathepsin L/V/K/H Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	Cathepsin L/V/K/H Antibody detects endogenous levels of total Cathepsin L/V/K/H
Immunogen Description	A synthesized peptide derived from human Cathepsin L/V/K/H
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Cathepsin; CTSL; CTSV; CTSK; CTSH;
Accession No.	Uniprot:P07711/O60911/P43235/P09668
Calculated MW	38kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Closely involved in osteoclastic bone resorption and may participate partially in the disorder of bone remodeling. Displays potent endoprotease activity against fibrinogen at acid pH. May play an important role in extracellular matrix degradation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.