

Dopamine Transporter Conjugated Antibody

Catalog No: #C56148

Package Size: #C56148-AF350 100ul #C56148-AF405 100ul #C56148-AF488 100ul #C56148-AF555 100ul #C56148-AF594 100ul #C56148-AF647 100ul #C56148-AF680 100ul #C56148-AF750 100ul #C56148-Biotin 100ul #C56148-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Dopamine Transporter Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	Dopamine Transporter Antibody detects endogenous levels of total Dopamine Transporter
Immunogen Description	A synthesized peptide derived from human Dopamine Transporter
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	DA transporter; DAT1; Dopamine transporter 1; PKDYS; SLC6A3;
Accession No.	Uniprot:Q01959
Calculated MW	50kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Amine transporter. Terminates the action of dopamine by its high affinity sodium-dependent reuptake into presynaptic terminals.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.