

## GRP94 Conjugated Antibody

Catalog No: #C56153



Package Size: #C56153-AF350 100ul #C56153-AF405 100ul #C56153-AF488 100ul #C56153-AF555 100ul #C56153-AF594 100ul #C56153-AF647 100ul #C56153-AF680 100ul #C56153-AF750 100ul #C56153-Biotin 100ul #C56153-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GRP94 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse
Specificity	GRP94 Antibody detects endogenous levels of total GRP94
Immunogen Description	A synthesized peptide derived from human GRP94
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Endoplasmin; GRP-94; gp96 homolog; HSP90B1; GRP 94; TRA1;
Accession No.	Uniprot:P14625
Calculated MW	94kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Molecular chaperone that functions in the processing and transport of secreted proteins. Functions in endoplasmic reticulum associated degradation (ERAD). Has ATPase activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.