

# Lactoferrin Conjugated Antibody

Catalog No: #C56155

Package Size: #C56155-AF350 100ul #C56155-AF405 100ul #C56155-AF488 100ul #C56155-AF555 100ul #C56155-AF594 100ul #C56155-AF647 100ul #C56155-AF680 100ul #C56155-AF750 100ul #C56155-Biotin 100ul #C56155-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Lactoferrin Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Lactoferrin Antibody detects endogenous levels of total Lactoferrin
Immunogen Description	A synthesized peptide derived from human Lactoferrin
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	GIG12; Growth inhibiting protein 12; HLF2; Kaliocin 1; Lactoferrin; Lactoferroxin-C; Lactotransferrin; LF; LTF; Neutrophil lactoferrin; Talalactoferrin;
Accession No.	Uniprot:P02788
Calculated MW	78kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Lactoferroxins A, B and C have opioid antagonist activity. Lactoferroxin A shows preference for mu-receptors, while lactoferroxin B and C have somewhat higher degrees of preference for kappa-receptors than for mu-receptors.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.