

Calpain 1 Conjugated Antibody

Catalog No: #C56175

Package Size: #C56175-AF350 100ul #C56175-AF405 100ul #C56175-AF488 100ul #C56175-AF555 100ul #C56175-AF594 100ul #C56175-AF647 100ul #C56175-AF680 100ul #C56175-AF750 100ul #C56175-Biotin 100ul #C56175-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Calpain 1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	Calpain 1 Antibody detects endogenous levels of total Calpain 1
Immunogen Description	A synthesized peptide derived from human Calpain 1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Calpain-1 catalytic subunit; CANP; CANP small subunit; CANPL1; CAPN1; Micromolar Calpain; Mu Calpain; muCANP;
Accession No.	Uniprot:P07384
Calculated MW	75,80kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

FC: 1:50-1:200

Product Description

Calcium-regulated non-lysosomal thiol-protease which catalyze limited proteolysis of substrates involved in cytoskeletal remodeling and signal transduction.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.