

# PAI1 Conjugated Antibody

Catalog No: #C56181

Package Size: #C56181-AF350 100ul #C56181-AF405 100ul #C56181-AF488 100ul #C56181-AF555 100ul #C56181-AF594 100ul #C56181-AF647 100ul #C56181-AF680 100ul #C56181-AF750 100ul #C56181-Biotin 100ul #C56181-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PAI1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	PAI1 Antibody detects endogenous levels of total PAI1
Immunogen Description	A synthesized peptide derived from human PAI1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	PAI-1; PAI-1; PAI; Endothelial plasminogen activator inhibitor; SERPINE1; PAI1; PLANH1; Serpine 1;
Accession No.	Uniprot:P05121
Calculated MW	45kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

This inhibitor acts as 'bait' for tissue plasminogen activator, urokinase, and protein C. Its rapid interaction with TPA may function as a major control point in the regulation of fibrinolysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.