

BNIP1 Conjugated Antibody

Catalog No: #C56182

Package Size: #C56182-AF350 100ul #C56182-AF405 100ul #C56182-AF488 100ul #C56182-AF555 100ul #C56182-AF594 100ul #C56182-AF647 100ul #C56182-AF680 100ul #C56182-AF750 100ul #C56182-Biotin 100ul #C56182-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	BNIP1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	BNIP1 Antibody detects endogenous levels of total BNIP1
Immunogen Description	A synthesized peptide derived from human BNIP1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	BCL2 adenovirus E1B 19kD interacting protein 1; BNIP1; Nip1; sec20; SEC20L; Transformation-related gene 8 protein; TRG8; Vesicle transport protein SEC20;
Accession No.	Uniprot:Q12981
Calculated MW	26kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

Product Description

SNARE that may be involved in targeting and fusion of Golgi-derived retrograde transport vesicles with the ER. Required for maintenance of ER network. Implicated in the suppression of cell death.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.