

# IRAK Conjugated Antibody

Catalog No: #C56188

Package Size: #C56188-AF350 100ul #C56188-AF405 100ul #C56188-AF488 100ul #C56188-AF555 100ul #C56188-AF594 100ul #C56188-AF647 100ul #C56188-AF680 100ul #C56188-AF750 100ul #C56188-Biotin 100ul #C56188-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	IRAK Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	IRAK Antibody detects endogenous levels of total IRAK
Immunogen Description	A synthesized peptide derived from human IRAK
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Il1rak; IRAK; Irak1; IRAK1-S; mPLK; Pelle; Pelle homolog; Pelle-like protein kinase; Plpk;
Accession No.	Uniprot:P51617
Calculated MW	77kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Binds to the IL-1 type I receptor following IL-1 engagement, triggering intracellular signaling cascades leading to transcriptional up-regulation and mRNA stabilization.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.