

Glypican 3 Conjugated Antibody

Catalog No: #C56189

Package Size: #C56189-AF350 100ul #C56189-AF405 100ul #C56189-AF488 100ul #C56189-AF555 100ul #C56189-AF594 100ul #C56189-AF647 100ul #C56189-AF680 100ul #C56189-AF750 100ul #C56189-Biotin 100ul #C56189-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Glypican 3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	Glypican 3 Antibody detects endogenous levels of total Glypican 3
Immunogen Description	A synthesized peptide derived from human Glypican 3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	DGSX; Glypican proteoglycan 3; Gpc3; GTR2 2; MXR7; OC15; SDYS; Secreted glypican-3; SGB; SGBS1;
Accession No.	Uniprot:P51654
Calculated MW	66kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Cell surface proteoglycan that bears heparan sulfate. Inhibits the dipeptidyl peptidase activity of DPP4. May be involved in the suppression/modulation of growth in the predominantly mesodermal tissues and organs.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.