

MTA2 Conjugated Antibody

Catalog No: #C56194

Package Size: #C56194-AF350 100ul #C56194-AF405 100ul #C56194-AF488 100ul #C56194-AF555 100ul #C56194-AF594 100ul
 #C56194-AF647 100ul #C56194-AF680 100ul #C56194-AF750 100ul #C56194-Biotin 100ul #C56194-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MTA2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	MTA2 Antibody detects endogenous levels of total MTA2
Immunogen Description	A synthesized peptide derived from human MTA2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Mata111; Metastasis associated gene 1 like 1; Metastasis associated gene family member 2; Mmta2; MTA1L1; p53 target protein in deacetylase complex; PID;
Accession No.	Uniprot:O94776
Calculated MW	75kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

FC: 1:50-1:200

Product Description

The p53 tumor suppressor gene integrates numerous signals that control cell life and death. There are several proteins that are involved in the p53 pathway, including Chk2, p53R2, p53AIP1, Noxa, PIDD, and MTA2. The transcriptional activity of p53 is modulated by protein stability and acetylation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.