

# CRMP2 Conjugated Antibody

Catalog No: #C56227

Package Size: #C56227-AF350 100ul #C56227-AF405 100ul #C56227-AF488 100ul #C56227-AF555 100ul #C56227-AF594 100ul #C56227-AF647 100ul #C56227-AF680 100ul #C56227-AF750 100ul #C56227-Biotin 100ul #C56227-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CRMP2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	CRMP2 Antibody detects endogenous levels of total CRMP2
Immunogen Description	A synthesized peptide derived from human CRMP2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	DRP2; CRMP-2; N2A3; Unc-33-like phosphoprotein 2; ULIP-2; DPYSL2; CRMP2; ULIP2; Musunc33; TOAD64;
Accession No.	Uniprot:Q16555
Calculated MW	62kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF: 1:50-1:200

FC: 1:50-1:200

## Product Description

Necessary for signaling by class 3 semaphorins and subsequent remodeling of the cytoskeleton. Plays a role in axon guidance, neuronal growth cone collapse and cell migration.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.