

Caspr Conjugated Antibody

Catalog No: #C56228

Package Size: #C56228-AF350 100ul #C56228-AF405 100ul #C56228-AF488 100ul #C56228-AF555 100ul #C56228-AF594 100ul #C56228-AF647 100ul #C56228-AF680 100ul #C56228-AF750 100ul #C56228-Biotin 100ul #C56228-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Caspr Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Caspr Antibody detects endogenous levels of total Caspr
Immunogen Description	A synthesized peptide derived from human Caspr
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Caspr; Caspr1; CNTNAP; Cntnap1; Contactin associated protein 1; MHDNIV; NCP1; Neurexin-4; Nrnx4; p190; Paranodin;
Accession No.	Uniprot:P78357
Calculated MW	165kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Contactin-associated protein 2 (Caspr2) is a type I transmembrane protein and member of the neurexin superfamily that mediates nervous system cell-cell interactions through the Neurexin IV-Caspr-Paranodin (NCP) complex.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.