

GCSF Receptor (CD114) Conjugated Antibody

Catalog No: #C56229



Package Size: #C56229-AF350 100ul #C56229-AF405 100ul #C56229-AF488 100ul #C56229-AF555 100ul #C56229-AF594 100ul #C56229-AF647 100ul #C56229-AF680 100ul #C56229-AF750 100ul #C56229-Biotin 100ul #C56229-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	GCSF Receptor (CD114) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	GCSF Receptor (CD114) Antibody detects endogenous levels of total GCSF Receptor (CD114)
Immunogen Description	A synthesized peptide derived from human GCSF Receptor (CD114)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CD114 antigen; Colony stimulating factor 3 receptor (granulocyte); CSF3R; Csfr; G-CSF receptor; GCSFR;
Accession No.	Uniprot:Q99062
Calculated MW	92kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Receptor for granulocyte colony-stimulating factor (CSF3), essential for granulocytic maturation. Plays a crucial role in the proliferation, differentiation and survival of cells along the neutrophilic lineage. In addition it may function in some adhesion or recognition events at the cell surface.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.