

Parathyroid Hormone Conjugated Antibody

Catalog No: #C56230

Package Size: #C56230-AF350 100ul #C56230-AF405 100ul #C56230-AF488 100ul #C56230-AF555 100ul #C56230-AF594 100ul #C56230-AF647 100ul #C56230-AF680 100ul #C56230-AF750 100ul #C56230-Biotin 100ul #C56230-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Parathyroid Hormone Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Parathyroid Hormone Antibody detects endogenous levels of total Parathyroid Hormone
Immunogen Description	A synthesized peptide derived from human Parathyroid Hormone
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	hPTH; Parathormone; Parathyrin; Parathyroid hormone 1; Prepro PTH; Preproparathyroid hormone; PTH1; PTH1 receptor; PTH1R; PTHR1;
Accession No.	Uniprot:P01270
Calculated MW	36kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

Product Description

PTH elevates calcium level by dissolving the salts in bone and preventing their renal excretion. Stimulates [1-14C]-2-deoxy-D-glucose (2DG) transport and glycogen synthesis in osteoblastic cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.