

Ran Conjugated Antibody

Catalog No: #C56233



Package Size: #C56233-AF350 100ul #C56233-AF405 100ul #C56233-AF488 100ul #C56233-AF555 100ul #C56233-AF594 100ul #C56233-AF647 100ul #C56233-AF680 100ul #C56233-AF750 100ul #C56233-Biotin 100ul #C56233-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Ran Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	Ran Antibody detects endogenous levels of total Ran
Immunogen Description	A synthesized peptide derived from human Ran
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ARA24; Gsp1; GTPase Ran; Guanosine triphosphatase Ran; LPS; RAN member RAS oncogene family; RanGTPase; Ras like protein TC4; Ras related nuclear protein; Ras-like protein TC4; RASL2 8;
Accession No.	Uniprot:P62826
Calculated MW	24kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

FC: 1:50-1:200

Product Description

GTP-binding protein involved in nucleocytoplasmic transport. Required for the import of protein into the nucleus and also for RNA export. Involved in chromatin condensation and control of cell cycle (By similarity).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.