

Emerin Conjugated Antibody

Catalog No: #C56235

Package Size: #C56235-AF350 100ul #C56235-AF405 100ul #C56235-AF488 100ul #C56235-AF555 100ul #C56235-AF594 100ul #C56235-AF647 100ul #C56235-AF680 100ul #C56235-AF750 100ul #C56235-Biotin 100ul #C56235-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Emerin Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Emerin Antibody detects endogenous levels of total Emerin
Immunogen Description	A synthesized peptide derived from human Emerin
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	EDMD; Emd; Emerin; Emery Dreifuss muscular dystrophy; STA;
Accession No.	Uniprot:P50402
Calculated MW	29kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Emerin may regulate gene expression through binding to other transcriptional regulators. Emerin binds to β -catenin and inhibits its nuclear accumulation. Recent studies demonstrate that Emerin is required for HIV-1 infectivity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.