

Dynamin 1 Conjugated Antibody

Catalog No: #C56249

Package Size: #C56249-AF350 100ul #C56249-AF405 100ul #C56249-AF488 100ul #C56249-AF555 100ul #C56249-AF594 100ul #C56249-AF647 100ul #C56249-AF680 100ul #C56249-AF750 100ul #C56249-Biotin 100ul #C56249-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Dynamin 1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Dynamin 1 Antibody detects endogenous levels of total Dynamin 1
Immunogen Description	A synthesized peptide derived from human Dynamin 1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	B dynamin; D100; DNM 1; DNM1; Dynamin; Dynamin1;
Accession No.	Uniprot:Q05193
Calculated MW	97kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Microtubule-associated force-producing protein involved in producing microtubule bundles and able to bind and hydrolyze GTP. Most probably involved in vesicular trafficking processes. Involved in receptor-mediated endocytosis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.