

Trefoil Factor 3 Conjugated Antibody

Catalog No: #C56250

Package Size: #C56250-AF350 100ul #C56250-AF405 100ul #C56250-AF488 100ul #C56250-AF555 100ul #C56250-AF594 100ul #C56250-AF647 100ul #C56250-AF680 100ul #C56250-AF750 100ul #C56250-Biotin 100ul #C56250-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Trefoil Factor 3 Conjugated Antibody |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB, IF |
| Species Reactivity | Human |
| Specificity | Trefoil Factor 3 Antibody detects endogenous levels of total Trefoil Factor 3 |
| Immunogen Description | A synthesized peptide derived from human Trefoil Factor 3 |
| Conjugates | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750 |
| Other Names | hITF; hP1.B; Intestinal trefoil factor; ITF; mITF; P1B; Polypeptide P1.B; TFF3; TFI; TREFOIL; Trefoil factor 3; |
| Accession No. | Uniprot:Q07654 |
| Calculated MW | 9kDa |
| Storage | Store at 4°C in dark for 6 months |

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Involved in the maintenance and repair of the intestinal mucosa. Promotes the mobility of epithelial cells in healing processes (motogen).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.