

CD5 Conjugated Antibody

Catalog No: #C56254

Package Size: #C56254-AF350 100ul #C56254-AF405 100ul #C56254-AF488 100ul #C56254-AF555 100ul #C56254-AF594 100ul #C56254-AF647 100ul #C56254-AF680 100ul #C56254-AF750 100ul #C56254-Biotin 100ul #C56254-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CD5 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	CD5 Antibody detects endogenous levels of total CD5
Immunogen Description	A synthesized peptide derived from human CD5
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CD5; CD5 antigen (p56 62); CD5 molecule; LEU1; LyA; Lymphocyte antigen T1/Leu-1; T-cell surface glycoprotein CD5; T1;
Accession No.	Uniprot:P06127
Calculated MW	55kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

Product Description

May act as a receptor in regulating T-cell proliferation. CD5 interacts with CD72/LYB-2.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.