

## CD82 Conjugated Antibody

Catalog No: #C56256



Package Size: #C56256-AF350 100ul #C56256-AF405 100ul #C56256-AF488 100ul #C56256-AF555 100ul #C56256-AF594 100ul #C56256-AF647 100ul #C56256-AF680 100ul #C56256-AF750 100ul #C56256-Biotin 100ul #C56256-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD82 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Rat
Specificity	CD82 Antibody detects endogenous levels of total CD82
Immunogen Description	A synthesized peptide derived from human CD82
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	C33; CD82; IA4; Inducible membrane protein; KAI1; Kangai1; SAR2; ST6; Tetraspanin 27; Tspan27;
Accession No.	Uniprot:P27701
Calculated MW	30-60kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Associates with CD4 or CD8 and delivers costimulatory signals for the TCR/CD3 pathway.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.