

MUM1 Conjugated Antibody

Catalog No: #C56258

Package Size: #C56258-AF350 100ul #C56258-AF405 100ul #C56258-AF488 100ul #C56258-AF555 100ul #C56258-AF594 100ul #C56258-AF647 100ul #C56258-AF680 100ul #C56258-AF750 100ul #C56258-Biotin 100ul #C56258-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	MUM1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	MUM1 Antibody detects endogenous levels of total MUM1
Immunogen Description	A synthesized peptide derived from human MUM1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	LSIRF; Multiple myeloma oncogene 1; NF-EM5; IRF4; MUM1; Transcriptional activator PIP;
Accession No.	Uniprot:Q15306
Calculated MW	52kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Melanoma associated antigen (mutated) 1 (MUM1, EXPAND1) is a PWWP-domain containing chromatin binding protein involved in maintaining chromatin architecture of interphase chromosomes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.