

KAT9 / Elp3 Conjugated Antibody

Catalog No: #C56261

Package Size: #C56261-AF350 100ul #C56261-AF405 100ul #C56261-AF488 100ul #C56261-AF555 100ul #C56261-AF594 100ul #C56261-AF647 100ul #C56261-AF680 100ul #C56261-AF750 100ul #C56261-Biotin 100ul #C56261-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	KAT9 / Elp3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	KAT9 / Elp3 Antibody detects endogenous levels of total KAT9 / Elp3
Immunogen Description	A synthesized peptide derived from human KAT9 / Elp3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Elongation protein 3 homolog; elp3; hELP3; Kat9;
Accession No.	Uniprot:Q9H9T3
Calculated MW	62kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Catalytic histone acetyltransferase subunit of the RNA polymerase II elongator complex, which is a component of the RNA polymerase II (Pol II) holoenzyme and is involved in transcriptional elongation. Elongator may play a role in chromatin remodeling and is involved in acetylation of histones H3 and probably H4. May also have a methyltransferase activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.