

# Otx1/2 Conjugated Antibody

Catalog No: #C56267

Package Size: #C56267-AF350 100ul #C56267-AF405 100ul #C56267-AF488 100ul #C56267-AF555 100ul #C56267-AF594 100ul #C56267-AF647 100ul #C56267-AF680 100ul #C56267-AF750 100ul #C56267-Biotin 100ul #C56267-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Otx1/2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Otx1/2 Antibody detects endogenous levels of total Otx1/2
Immunogen Description	A synthesized peptide derived from human Otx1/2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Homeobox protein OTX1; Homeobox protein OTX2; Orthodenticle 1; Orthodenticle 2; Orthodenticle homeobox 1; Orthodenticle homeobox 2; Orthodenticle1; Orthodenticle2; otx1; otx2;
Accession No.	Uniprot:P32242,P32243
Calculated MW	32kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Probably plays a role in the development of the brain and the sense organs. Can bind to the BCD target sequence (BTS): 5'-TCTAATCCC-3'.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.